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Saw blades, non sterile – Instructions for use and cleaning

EN

Products
Saw blades for bone surgery, non-sterile

REF number groups

10xxx	11xxx	12xxx	13xxx	14xxx	15xxx
16xxx	17xxx	18xxx	19xxx	20xxx	21xxx
22xxx	23xxx	24xxx			

By purchasing this instrument, you have chosen a high-quality product whose proper use and handling are illustrated in the following. To keep the risks for patients and users to a minimum, we ask you to read and follow the instructions for use carefully. The instruments may only be used, disinfected, cleaned and sterilized by trained and qualified professionals.

Purpose The products are designed to cut bones in combination with an electric drive.

Qualification of the user These instruments are only designed for qualified doctors with sufficient experience in the field of bone surgery. Explanations of detailed surgical procedures are therefore omitted.

Tests / function control The instruments must be identified and checked for proper function before each use. The instruments may only be used sterilized. Damage to the surfaces such as scratches, cracks, nicks, notches etc., as well as bent parts mean that the product must not be used. In this case, the products must be disposed of in the normal manner for hospitals.

Do not use damaged products!
(see tests / function control)

Defective products must be subjected to the full reprocessing procedure before they are sent in for repair. Products that cannot be repaired any more are to be disposed of in the normal manner for hospitals.

Compatibility Depending on the connection, the saw blades can be operated with the appropriate drive machines or commercially available motor systems. The compatibility does not depend on the machine type but on the connections of the saw blade. The connections of the saw blades listed below are compatible with the respective devices of the manufacturers mentioned: (risa GmbH is not affiliated with these companies.)

REF (number group)	Description / Power Tool
10xxx	Saw blade, Synthes / AO-coupling
11xxx	Saw blade, Aesculap-coupling
12xxx	Saw blade, Linvatec-coupling
13xxx	Saw blade, Linvatec / 3M-coupling
14xxx	Saw blade, Stryker-coupling
15xxx	Saw blade, Micro-Aire-coupling
16xxx	Saw blade, DeSoutter-coupling
17xxx	Saw blade, Tava-coupling
18xxx	Saw blade, Dyonics-coupling / (Tuke-Saw)
19xxx	Saw blade, Dreikant-coupling
20xxx	Saw blade, Arthrex-coupling
21xxx	Saw blade, Bojin-coupling
22xxx	Saw blade, custom made
23xxx	Saw blade, Medicon / Nouvag-coupling
24xxx	Saw blade, Arthrex-coupling

Area of application
Indication Saw blades to cut bones.
Contraindication The products are contraindicated for all applications beyond the intended purpose.
Patient group There are no restrictions on age or aetiology.
User group The products are only designed for qualified doctors with sufficient experience in the field of bone surgery. Explanations of detailed surgical procedures are therefore omitted.
Mode of operation The saw blades are connected to a power tool (sawing machine) via the connection at the proximal end. Activating this drive unit causes the saw blade to oscillate in a controlled manner. The saw blade, which has sharp teeth, gradually removes bone material through the oscillation until the desired depth of cut is reached or the bone is severed.

Handling
Do not use damaged products!
(see tests / function control)

The saw blades may only be used with a power tool that matches the connection, see compatibility.
Before use, ensure that it is seated correctly in the holder (the tool must not have any play in the holder). The saw blade must be driven before attaching to the bone. The lowest possible pressure applied by the saw blade, a high oscillation speed and the use of instruments that are as sharp as possible reduce the risk of overheating and the associated thermal necrosis. Avoid contact with metallic objects as the saw blades can become blunt. The saw blades must not be overstressed by twisting or prying, as this can damage or break instrument parts. The saw blades contain no user-serviceable parts. The saw blades must not be repaired – this can lead to injuries to the user or patient.

Risks
- Mechanical risks (break, fragmentation) due to misuse (e.g. too high pressing force) or wearing
- Thermic risks (tissue necrosis due to temperature increase in surrounding tissue), misuse or wearing
- Biological / microbiological risks due to inaccurate preparation by the user

Disposal If the instruments cannot be repaired and reprocessed any more, then they should be disposed of in the normal manner for hospitals.

Materials The materials used are stainless steels according to DIN EN ISO 7153-1.

Instructions for reprocessing according to DIN EN ISO 17664-1 Method: Cleaning, Disinfection, Sterilization hot steam (DIN EN ISO 17665-1)

Notes The instruments are delivered non-sterile and must be cleaned and possibly disinfected and sterilized before use. A cleaning - / sterilization method was validated according to DIN EN ISO 17664 and will be illustrated in the following. The efficiency and suitability of the procedures below for the preparation of the products concerned has been validated by the manufacturer. It is the responsibility of the operator to ensure effective preparation using appropriate equipment, process aids and adequately qualified personnel for implementation. Any deviation from the procedures specified below must be checked by the user / operator for its influence regarding the safe and effective preparation of the products. In the event of deviations from the procedures described below, the manufacturer assumes no responsibility for the effectiveness and compatibility of these procedures. The saw blades may only be prepared by persons who have the necessary specialist knowledge and training and who can assess the risks involved with the corresponding effects.

Restriction of reprocessing Due to the design of the products and the materials used, we cannot specify a defined limit in terms of the maximum number of reprocessing cycles. The service life of the medical products is determined by their function and how carefully they are handled and used.

Do not use any damaged products!
(see tests / function control)

Place of use In order to prevent incrustation of dirt and coagulation of blood residues, visible dirt on the product should be removed immediately after use. Avoid temperatures above 40 °C, as these can promote the fixation of contamination. The products are to be transferred in a suitable, closed container within one hour of processing. Immediately after use, instruments should be manually disinfected in order to reduce the risk of an infection for the user.
To disinfect the instruments, place them in disinfectant solution. It must be ensured that the instruments are completely submerged in the disinfectant solution and that no air bubbles have formed on them. Follow the instructions of the manufacturer of the disinfectant solution.

Preparation for decontamination If the instruments can be taken apart, then take them apart before reprocessing.

Manual pre-cleaning
1. Rough Pollution Remove rough dirt under cold tap water with an appropriate brush or cloth.
2. Cleaning bath Place the products into a solution with 0.5 – 2% Neodisher® MediClean forte with tap water (< 30 °C, drinking water quality) for 10-30 minutes. The surface and cavity must be fully wetted. Then clean the surfaces, channels, lumens and cavities of the instrument under running cold water (< 30 °C, drinking water) using a suitable brush (no steel brush). Flush through gaps, slits and cavities for about 10 seconds with a water pressure gun.
3. Ultrasonic bath Place instruments in an appropriate ultrasonic bath. Instruments must be fully wetted with cleaning solution (e.g. 0.5 - 2% Neodisher® mediclean forte in demineralized water). Beware of using appropriate sterilization trays. Do not overload the trays. Make sure that there are no sonic shadows. The ultrasonic bath must be replaced before each use. Clean parts for 10 minutes in ultrasonic bath at room temperature and at a frequency of a minimum of 35 kilohertz.
4. Rinse (This step can be omitted if automatic cleaning is immediately continued.)
 Rinsing the entire instrument with demineralized water for at least 1 minute.
5. Drying (This step doesn't apply if cleaning is immediately continued.)
 Each instrument has to be entirely dry inside and outside after pre-cleaning to prevent rust creation or mal-function. Dry the instruments with a soft cloth / swab. A compressed air system for drying can be used to dry hollow parts or the inner channel.

Automatic cleaning and disinfection
A mechanical method (cleaning and disinfection device) should be used to clean / disinfect the instruments.
The following points in particular must be noted for mechanical reprocessing:
- A prerequisite for effective mechanical reprocessing is the proper loading of perforated instrument trays, inserts, holders etc. The perforated instrument trays must not be overloaded so that the instruments can be washed well.
- Large instruments must be placed on the perforated instrument tray so that they do not prevent other instruments from becoming clean.
- Instruments with cavities must also be completely washed on the inside.
- The instruments must be placed according to their mechanical sensitivity so that damage is excluded. Sharp cutting edges must be exposed.
- The washed instruments must be removed immediately from the machine after the cleaning program finishes because leaving them in the closed machine can cause corrosion due to moisture.

Optional:
Thermal disinfection after machine cleaning: min. 5 min at 92 °C (+/- 2 °C).
Thermal disinfection per se is not a final process that puts the product in a ready-to-use condition. This disinfection step only improves infection protection between cleaning and sterilization. A subsequent sterilization according to the procedure described below is always necessary.

Recommended cleaning method
A mechanical reprocessing program with thermal disinfection includes the following steps:
Cleaning and disinfection unit: Miele G 7835 CD
Cleaning program: Vario TD AD
1. Pre-wash with cold tap water (< 30 °C, drinking water quality)
2. Drain water
3. Pre-wash with cold tap water
4. Drain water
5. Cleaning with 0,5 % Neodisher® Mediclean forte in deionised water at a temperature of 55 °C (+/- 1 °C) for minimum 10 minutes.
6. Drain water
7. Rinse 3 minutes with cold tap water
8. Drain water
9. Rinse 2 minutes with deionised water
10. Thermal Disinfection at 92 ° C (+/- 2 °C) > 5 minutes
11. Dry the instruments in the drying cycle of the cleaning / disinfection device at 100 °C for 30 minutes. If necessary, manual drying can also be performed with the help of a lint-free cloth. Dry cavities in instruments with sterile compressed air.
 Inadequately cleaned instruments must be cleaned, disinfected and adequately rinsed off again as described under „Automatic cleaning and disinfection". According to the RKI guideline "Hygiene Requirements for the Reprocessing of Medical Devices", a manual method is not specified and validated.

Packaging To comply with standards, use (double) packaging for the instruments for sterilization in sterile bags conforming to ISO 11607-1. Deviations from the suggested packaging must be validated.

Sterilization Recommended sterilization method: Sterilization parameters (steam sterilization with moist heat according to DIN IN ISO 17665-1, type B autoclave):
at least triple pre-vacuum, 134 °C for 5 minutes, one drying cycle for at least 20 minutes. Check the sterile goods packaging for damage after sterilization and check the sterilization indicators.

Storage Store sterilized instruments in sterile goods packing in a closed cabinet protected from dust, moisture and temperature fluctuations.

Additional information If the chemicals and machines described here are not available and if the reprocessing procedure cannot be performed as described, then it is the responsibility of the user to validate the corresponding method used.
Additional information on the reprocessing of medical devices:
- Internet: http://www.rki.de
- Internet: http://www.a-k-i.org
- Recommendations of the Commission for Hospital Hygiene and Infectious Desease Prevention of the Robert Koch Institute (RKI) and of the German Federal Institute for Drugs and Medical Devices (BfArM) regarding the "Hygiene requirements for the Reprocessing of Medical devices".

Warranty The products are manufactured from high quality materials and are subject to quality control before delivery. However, if defects are found in spite of this, then please contact our service department. risa GmbH cannot provide any guarantee that the products are suitable for the particular surgical procedure. The users must decide this for themselves.
risa GmbH is not liable for accidental or consequential damage. risa GmbH cannot be held liable in cases in which it is proven that the instructions for use were not followed properly.
Caution:
In case where the instruments are used on patients with Creutzfeldt-Jakob disease or an HIV infection, we dis-claim any responsibility for reprocessing.
To process the instruments after treating cases of Creutzfeldt-Jakob disease or an HIV infection, special validations and routine monitoring of the method by the user are necessary.

risa GmbH
Industriestr. 7
78234 Engen, Germany
info@risa-germany.eu
www.risa-germany.eu

CE 0483

Symbols

	Manufacturer	REF	Article number
LOT	Lot number-		Follow instructions for use
NON STERILE	non sterile		Caution